SEALTEST

PT

TESTE DE SELAGEM PARA SELADORAS

SEALTEST - EMBALAGEM COM 100 UNIDADES

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EN SEALING TEST FOR SEALER

PRODUCT DESCRIPTION

The **SEALTEST MAXXIMED** evaluates the correct functioning of automatic sealers in terms of the quality of the heat sealing performed when closing medical packages made of surgical grade paper that contain health products.

Sealers must promote heat-sealing without cracks, wrinkles or signs of burning of the packages, because these occurrences can cause difficulties in the aseptic transfer of the content or delaminations that release particles when the package is opened and can contaminate the health product or the surgical field.

It must be part of the CME routine as a vital instrument to ensure the hermetic sealing of packages, resulting in a sterile barrier system that allows sterilization, provides an appropriate microbial barrier and allows the aseptic presentation of the contents from sterilization, through storage and transport, until the moment of use.

The **SEALTEST MAXXIMED** is available in boxes of 100 units, in accordance with the ABNT NBR ISO 11607-2 standard.

COMPOSITION

Each **SEALTEST MAXXIMED** unit is produced on a sheet of paper printed with thermo-reactive ink that evaluates the sealing parameters (temperature, pressure and sealing speed). It makes it possible to identify abnormalities in the operation of sealers and prevent packaging losses due to inefficient sealing, generating savings and productivity.

In dimensions 20 cm x 6 cm (W x H), the product is entirely made of recyclable materials, endorsing the MAXXIMED commitment to the conscious use of resources for the preservation of the environment.

INSTRUCTIONS FOR USE

1. Turn on and wait for the sealer to reach the ideal temperature of operation.
2. Fill in the data contained in the **SEALTEST MAXXIMED** sheet.
3. Place the test sheet inside an empty container of sanitary grade paper, between the laminated plastic film and the paper, with the test side facing the film.
4. Carry out a sealing of the packaging, so that it horizontally covers the entire length of the test area.
5. At the end, go over the sealing aspects in detail and consider the following:
 - a) **Approved Result:** intact sealing of the packaging, without cracks, wrinkles or burning of the paper and test sheet without bubbles or stains.
 - b) **Disapproved Result:** detachment of the sealing on the packaging, presence of cracks, wrinkles or burning of the paper or test sheet with bubbles or stains.

IMPORTANT

- If approved, release the equipment for use.

- In the event of disapproval, repeat the test and if it results in approval, release it for use. If the failed/negative result repeats, interdict the equipment and call Maintenance / Engineering for the appropriate measures.

- Record the results and archive them for quality management and/or auditing purposes.

STORAGE

Store away from light in a clean, dry and ventilated place, at temperatures between 5°C and 40°C and relative humidity between 30% and 60%. Do not leave in direct contact with the floor, wall, chemical or sterilizing agents.

VALIDITY

The **SEALTEST MAXXIMED** is valid for 05 years from the date of manufacture. Lot information, expiration date and other identification data are printed on the packaging and on each unit of the product.

STANDARDS

The **SEALTEST MAXXIMED** is in accordance with ABNT NBR ISO 11607-2. Mandatory product to evaluate the operation of automatic sealers.

BIO PROGRAM

Green consciousness is the logic that drives MAXXIMED's business model. For this reason, the company is proud to have its people, facilities, processes, products and services incorporated into a sustainable system.

GREEN PARTNER PROGRAM

To ensure that the components used in its products are harmless to the environment, MAXXIMED incorporates the "Green Partner" into its processes, a program that approves the company's suppliers and partners according to appropriate environmental parameters and encourages the conduct of business that generates value over time.

LEAD FREE

MAXXIMED uses lead-free technologies and other heavy materials that dissolve when in contact with acidic water, contaminating water, air and soil, seriously harming human health.

FREE FROM VOLATILE ORGANIC COMPOUNDS

MAXXIMED uses technologies free from Volatile Organic Compounds that, when exposed to closed places, can worsen human health.

These are some of the actions in which MAXXIMED reinforces its commitment to continuous improvement and the preservation of the environment. Copyright 2019. Clicktest, Maxitest and Maxipack are registered trademarks of Maximed Latin America. Reproduction in whole or in part without written permission is prohibited.