
LIQUID 48H

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EN BIOLOGICAL INDICATORS FOR LIQUID LOAD STEAM 48H (IBVCL48H)

PRODUCT DESCRIPTION

The **LIQUID Colorimetric Biological Indicator – IBVCL48H (STEAM)** was developed to monitor steam sterilization cycles of liquid loads, between 121°C and 134°C. It makes it possible to efficiently prove the microbiological lethality of sterilization cycles.

The **LIQUID Colorimetric Biological Indicator – IBVCL48H (STEAM)** is available in boxes of 50 units and complies with AAMIISO 11138-1, AAMIISO 11138-3 and USP standards.

COMPOSITION

The **LIQUID Colorimetric Biological Indicator – IBVCL48H (STEAM)** is a self-contained biocidator made with a hermetically sealed glass ampoule, containing a minimum population of 1×10^6 of *Geobacillus stearothermophilus* (ATCC 7953) spores arranged in a purple-colored culture medium and pH indicator.

INSTRUCTIONS FOR USE

1. Remove the **LIQUID Colorimetric Biological Indicator – IBVCL48H (STEAM)** from the refrigerator and leave it at room temperature for 10 minutes.

2. Place the biological indicator inside a container with the liquid similar to the one that will be sterilized, in an area considered critical (of greatest challenge) of the sterilizer, as described in the sterilization validation protocol.

Observation: if there are containers and/or liquids and/or volumes of different configurations, other biological indicators must be used to monitor them.

3. Activate the sterilization cycle corresponding to the load to be processed. At the end of the cycle, open the sterilizer door, carefully remove the container with the biological indicator and wait until it cools down (approx. 5 minutes) and then, remove the biological indicator.

WARNING: use safety glasses and gloves when removing the **LIQUID Colorimetric Biological Indicator – IBVCL48H (STEAM)** from the sterilized container.

4. Take the processed biological indicator and a **LIQUID Colorimetric Biological Indicator – IBVCL48H (STEAM)** that has **NOT been subjected to the sterilization process** to be used as a positive control, in order to validate the incubation conditions (temperature and time), the viability of the spores and the qualification of the culture medium. Both biological indicators (processed and positive control) must belong to the same manufacturing lot.

5. Check that the incubator is operating at 60°C (+/- 2°C) and incubate the two biological indicators for 48 hours. Considering that the presence of microorganisms can be detected before 48 hours, it is recommended to take readings every 12 hours. Results:

a) maintenance of the purple color of the culture medium indicates the absence of viable microorganisms, therefore, a successful sterilization process and **satisfactory result**.

b) change in the color of the culture medium, from purple to yellow, indicates the presence/survival of microorganisms, therefore, sterilization cycle failure and **unsatisfactory results**.

6. Release the load for use only if the **processed LIQUID Colorimetric Biological Indicator – IBVCL48H (STEAM)** presents a result "Negative - Absence of Microorganisms" and the control indicator results "Positive - Presence of microorganisms", confirming that the sterilization cycle was satisfactory and that the incubation was adequate.

Observation: if the processed biological indicator shows a result "Positive - Presence of microorganisms" and, consequently, an irregularity in the sterilization cycle, the Hospital Infection Control Commission and Engineering/Maintenance must be notified to adopt the appropriate measures.

7. Record in a book the results obtained.

8. Dispose of the biological indicator ampoules following the sanitary regulations in your region. Positive indicators can be sterilized in cycles of 121°C for 15 minutes or 132°C for 10 minutes in gravitational Auto-Reader Incubators; or in cycles of 134°C for 04 minutes in vacuum pump-assisted Auto-Reader Incubators.

ATTENTION: do not use the **LIQUID Colorimetric Biological Indicator – IBVCL48H (STEAM)** to monitor ethylene oxide sterilization cycles, dry heat, plasma/hydrogen peroxide vaporization, low-temperature steam and formaldehyde, chemical vapors, radiation or other sterilization processes that do not use saturated steam as a sterilizing agent. Do not reuse biological indicators.

IMPORTANT: reading by colorimetric methodology in 48 hours meets the requirements of the FDA's reduced incubation time protocol (RIT), and the results are correlated with the incubation time of 7 days, in 97% to 100% of cases.

VALIDITY/ANVISA REGISTRATION

The **LIQUID Colorimetric Biological Indicator – IBVCL48H (STEAM)** has validity of 02 years from the date of manufacture. Information about the lot, date of manufacture and expiration date is printed on the packaging and on the biological indicator ampoule. Product not considered a healthcare product and exempt from registration with ANVISA.

STORAGE

Store away from light, in a refrigerator at temperatures between 2°C and 8°C and relative humidity between 30% and 80%. Do not freeze. Do not store in contact with floors and walls or near heat sources, sterilizing agents and other chemical products. Storage conditions must be strictly followed.

TRANSPORT

Transport in sealed and reinforced boxes to avoid damage. Transporting this product does not represent a risk to health.

ATTENTION: the ampoules can be kept at ambient temperatures below 25°C for a period of up to 15 days. If the product is exposed to temperatures **NOT recommended**, the stability of the contents of the ampoule may be compromised.