
# MAXITEST TIPO 6

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## EN CHEMICAL INDICATOR - TYPE 6 EMULATOR STEAM

### PRODUCT DESCRIPTION

The **MAXITEST CHEMICAL EMULATOR - Type 6** is a chemical indicator for internal use developed to carry out a complete and integrated assessment of the critical parameters of the sterilization cycle (steam quality, exposure time, and temperature). It promotes qualified reading of the conditions necessary for sterilization inside each package to be processed.

In compliance with the AAMI/ISO 11140-1 standard, the product guarantees efficiency and safety in the monitoring of steam sterilization cycles and its main feature is to simulate a specific cycle.

The **MAXITEST CHEMICAL EMULATOR - Type 6** is made of a medical grade paper strip, laminated front and back, and printed with thermochromic chemical formula reacting to saturated steam. Provides result reading by sharp color change that simplifies interpretation.

Dimensions: 19mm x 78mm (H x W). Presentation: boxes with 250 units.

### TECHNICAL INFORMATION

According to the definition and qualification requirements contained in the AAMI/ISO 11140-1 standard, Type 6 chemical indicators are called "cycle checkers" and are highly accurate in reading the declared values of temperature and time of exposure to steam in sterilization cycles. The **MAXITEST CHEMICAL EMULATOR - Type 6** fully meets these performance requirements and is calibrated to monitor the sterilization cycles of 121°C for 20 minutes and 134°C for 04 minutes.

This product is an essential element in the composition of the Quality Management System in the reprocessing of materials and medical devices.

### INSTRUCTIONS FOR USE

1. Insert one unit of **MAXITEST CHEMICAL EMULATOR - Type 6** inside each unit/package to be sterilized. Another **MAXITEST CHEMICAL EMULATOR - Type 6** must also be used in the assembly of a challenge test package that will work as a load releaser, according to ANVISA RDC 15: 2012.

2. Pack the load inside the auto-reading incubators' chamber, respecting the configuration approved in the validation (type of load, volume, and weight).

3. Select and activate the sterilization cycle corresponding to the load to be sterilized.

4. At the end of the cycle, carefully and respecting the auto-reading incubators' cooling time recommended by the manufacturer, remove the **MAXITEST CHEMICAL EMULATOR - Type 6** from the inside of the test package (if it has been used) and read the result.

Satisfactory result: and load released for storage in arsenal or for use in units: when the dot (circle) of the chemical indicator ink changes completely and uniformly from the original purple color to shades of brown color to black color. Compare the change to the Interpretation Guide printed on the strip and in this document.

**IMPORTANT:** consider as an accepted/approved result the uniform change from purple to brown or black shades. This situation occurs as a function of the conditions generated in the conditioning phase (vacuum) and/or the time spent in the heating ramp phase inside the autoclave chamber (time that takes from the end of the vacuum phase until reaching the temperature of sterilization). These phases expose the **MAXITEST CHEMICAL EMULATOR - Type 6** to conditions that exceed those established by the ISO 11140 standard and are possible only in resistors (Bier Vessel). Thus, variations in shades of brown and black are acceptable due to the non-uniformity of steam autoclave sterilization cycles and different load conditions.

5. If an unsatisfactory result is evidenced and the sterilization cycle fails, all materials contained in the load must be reprocessed.

6. The chemical indicator is a legal document. Therefore, its result must be recorded in an appropriate book/control and be attached to the patient's medical record, according to the routine established in the institution and health legislation in force.

**CAUTION:** Do not use the **MAXITEST CHEMICAL EMULATOR - Type 6** to monitor ethylene oxide sterilization cycles, dry heat, chemical vapors, radiation, or other sterilization processes that do not use saturated steam as a sterilizing agent.

### VALIDITY/ANVISA REGISTRATION

The **MAXITEST CHEMICAL EMULATOR - Type 6** is valid for 03 years from the date of manufacture. Lot information, expiration date and other identification data are printed on the packaging and on each unit of the product.

Product exempt from registration with ANVISA, according to RDC 185/2001.

### STORAGE

Store away from light at temperatures between 15°C and 30°C and relative humidity between 30% and 60%. Do not freeze. Do not store in contact with floors and walls or near heat sources, sterilizing agents and other chemical products.

### BIO PROGRAM

Green consciousness is the logic that drives MAXIMED's business model. For this reason, the company is proud to have its people, facilities, processes, products and services incorporated into a sustainable system.

### GREEN PARTNER PROGRAM

To ensure that the components used in its products are harmless to the environment, MAXIMED incorporates the "Green Partner" into its processes, a program that approves the company's suppliers and partners according to appropriate environmental parameters and encourages the conduct of business that generates value over time.

### LEAD FREE

MAXIMED uses lead-free technologies and other heavy materials that dissolve when in contact with acidic water, contaminating water, air and soil, seriously harming human health.

### FREE FROM VOLATILE ORGANIC COMPOUNDS

MAXIMED uses technologies free from Volatile Organic Compounds that, when exposed to closed places, can worsen human health.

These are some of the actions in which MAXIMED reinforces its commitment to continuous improvement and the preservation of the environment. Copyright 2019. Clicktest, Maxitest and Maxipack are registered trademarks of Maximed Latin America. Reproduction in whole or in part without written permission is prohibited.