
CLICKTEST 8H

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▶ REFERÊNCIAS

- Plug, J.J., Microbiology and Engineering of Sterilization Processes. 7th edition, Published by the Environmental Sterilization Laboratory, Minneapolis, MN, 1990.
- ANSI/AAM/ISO 11138-1:2006, Sterilization of health care products Biological Indicators Part 1: General requirements.
- EN 866-3:2000, Biological systems for testing sterilizers and sterilization processes Part 3: particular systems for use moist heat sterilizers.

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EN BIOLOGICAL INDICATORS - STEAM 8H (IBVC8H)

► PRODUCT DESCRIPTION

The **Clicktest Colorimetric Biological Indicator - IBVC8H STEAM** was developed to monitor sterilization cycles from 121°C to 134°C in vacuum pump assisted Auto-Reader Incubators, with air removal by gravity and immediate use cycle. It makes it possible to efficiently and safely prove the microbiological lethality of sterilization cycles.

Its innovative design dispenses glass components and, combined with an exclusive activation system, simplifies use and provides safe handling by preventing any risk of accidents.

The **Clicktest Colorimetric Biological Indicator - IBVC8H STEAM** is available in boxes of 50 units according to the AAMIISO 11138-1, AAMIISO 11138-3 and AAMIISO 11140-1 standards.

► COMPOSITION

Each **Clicktest Colorimetric Biological Indicator - IBVC8H STEAM** is made with thermoplastic tube, stainless steel ACS, special filter, adhesive label with process chemical indicator (Type 1), minimum population of 1×10^6 of *Geobacillus stearothermophilus* spores (ATCC 7953) inoculated in spore-bearing strip. It also has a purple-colored growth medium contained in a casing and thermoplastic lid with holes, plus a vapor permeable barrier.

► INSTRUCTIONS FOR USE

1. Identify the **Clicktest Colorimetric Biological Indicator - IBVC8H STEAM** by writing the sterilizer number, cycle/load number, and processing date on its label.

2. Pack the biological indicator in a standardized package according to recommended sterilization practices (ANMISA, SOBECC, APECIH, etc.) and institution validated.

3. Place the package next to the materials to be sterilized, in an area considered critical (of greatest challenge) of the Auto-Reader Incubators, usually located near the door and next to the drain.

4. Select and activate the sterilization cycle corresponding to the load to be processed. At the end of the cycle, wait for the cooling time recommended by the Auto-Reader Incubators manufacturer, open the chamber and carefully remove the biological indicator from the inside package.

5. Verify that the Type 1 external chemical indicator (contained on the label) has changed from its pink (original) color to brown or black, confirming that it has been exposed to vapor at 121°C or above.

Important: The color change is not evidence that the sterilization cycle was successful. If there is no color change, it will be necessary to review the sterilization process.

6. Activate the **Clicktest Colorimetric Biological Indicator - IBVC8H STEAM**. To do this, hold the ampoule and press the lid down **until you hear a "click"**, breaking the retaining seal of the culture medium. Shake it **until all the liquid** (culture medium) **has completely descended** from the lid and comes into contact with the spore-bearing strip (see Quick Manual in this document).

7. Take another **Clicktest Colorimetric Biological Indicator - IBVC8H STEAM** that **has NOT been subjected to the sterilization process** to be used as a positive control, in order to validate the incubation conditions (temperature and time), the viability of the spores and the qualification of the culture medium. Both biological indicators (processed and positive control) must belong to the same manufacturing lot. Repeat step 1 (identification) and step 6 (activation).

8. Check that the incubator is operating between 55°C and 60°C and incubate the 02 biological indicators for 08 hours to obtain the result by visual reading through colorimetric methodology (change of color of the culture medium).

Since the presence of microorganisms can be detected before 24 hours, it is recommended to take readings every 04 hours. Results:

a) maintenance of the purple color of the culture medium indicates the absence of viable microorganisms. Therefore, a successful sterilization process and **satisfactory result**.

b) change in the color of the culture medium, from purple to gray/greenish and later to yellow, indicates the presence/survival of microorganisms. Therefore, sterilization cycle failure and **unsatisfactory results**.

9. Release the load/load for use only if the processed **Clicktest Colorimetric Biological Indicator - IBVC8H STEAM** presents a result "Negative - Absence of Microorganisms" and the **control indicator** results "Positive - Presence of microorganisms", confirming that the sterilization cycle was satisfactory and that the incubation was adequate.

Note: if the processed biological indicator shows a result "Positive - Presence of microorganisms" and, consequently, an irregularity in the sterilization cycle, the Hospital Infection Control Commission and Engineering/Maintenance must be notified to adopt the appropriate measures.

10. Detach the biological indicator labels (processed and positive control), attach them to the designated book/file and record the results obtained.

11. Dispose of the biological indicator tubes/ampoules following the sanitary regulations in your region. Positive indicators can be sterilized in cycles of 121°C for 15 minutes or 132°C for 10 minutes in gravitational Auto-Reader Incubators; or in cycles of 132°C for 04 minutes in vacuum pump-assisted Auto-Reader Incubators.

All materials that make up the **Clicktest Colorimetric Biological Indicator - IBVC8H STEAM** are 100% recyclable, following the MAXXIMED institutional plan of ecological design.

CAUTION: Do not use the **Clicktest Colorimetric Biological Indicator - IBVC8H STEAM** to monitor ethylene oxide sterilization cycles, dry heat, plasma/hydrogen peroxide vaporization, low-temperature steam and formaldehyde, chemical vapors, radiation, or other sterilization processes that do not use saturated steam as a sterilizing agent. Do not reuse biological indicators.

► VALIDITY/ANVISA REGISTRATION

The **Clicktest Colorimetric Biological Indicator - IBVC8H STEAM** has validity of 02 years from the date of manufacture. Information about the lot and due date are printed on the packaging and on the biological indicator tube. Product exempt from registration with Anvisa, according to RDC 185/2001.

► STORAGE

Store away from light at temperatures between 15°C and 30°C and relative humidity between 35% and 60%. Do not freeze. Do not store in contact with floors and walls or near heat sources, sterilizing agents, and other chemical products.

IMPORTANT: always transport, store, activate, incubate and discard in an upright position (lid up) to prevent leakage after product activation.