
# MAXITEST™ TIPO 4

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## MULTIPARAMETRIC - HYDROGEN PEROXIDE (VH2O2)

### PRODUCT DESCRIPTION

The **MAXITEST CHEMICAL INDICATOR - Type 4 (VH2O2)** is an internal chemical indicator categorized as a multiparametric designed to evaluate plasma sterilization or hydrogen peroxide vaporization cycles. It promotes qualified reading of the conditions necessary for sterilization (hydrogen peroxide concentration, exposure time and temperature) not inside each package to be processed.

Produced in a controlled environment and in compliance with the ANSI/AAMIISO 11140-1 standard, the product guarantees efficiency and safety in monitoring sterilization cycles. Its result reading by color change simplifies and facilitates interpretation.

In dimensions 19mm x 78mm (H x W) in boxes with 250 units, the **MAXITEST CHEMICAL INDICATOR - Type 4 (VH2O2)** is made with a plastic strip, printed with an exclusive chemical formulation that reacts to the critical variables of the plasma or vaporization sterilization process, of hydrogen peroxide (VH2O2).

### TECHNICAL INFORMATION

The **MAXITEST CHEMICAL INDICATOR - Type 4 (VH2O2)** fully complies with the performance requirements contained in the ANSI/AAMIISO 11140-1 standard and is calibrated to react to the following parameters: concentration of the sterilizing agent (2,3mg/L), time exposure time (6 minutes) and exposure temperature (50°C)

It is an essential element in the composition of the Quality Management System and must be used in the routine of reprocessing medical materials and devices.

### INSTRUCTIONS FOR USE

1. Insert one unit of **MAXITEST CHEMICAL INDICATOR - Type 4 (VH2O2)** not inside each package that packs a product for health.

**NOTE:** to increase the quality of process monitoring, it is recommended to use another unit of the **MAXITEST CHEMICAL INDICATOR - Type 4 (VH2O2)** to assemble a challenge test package, which will work as a load release.

2. Pack the load and the challenge test pack inside the sterilizer chamber, respecting the equipment manufacturer's recommendations and the load configuration approved in validation (type of material, quantity, volume and weight).

3. Select and activate the sterilization cycle corresponding to the load to be sterilized.

4. At the end of the cycle, carefully and respecting the manufacturer's recommendations, remove the **MAXITEST CHEMICAL INDICATOR - Type 4 (VH2O2)** from the inside of the challenge test package (if it has been used) and read the result.

Satisfactory result and load released for storage in the armory or for use in units: when the chemical indicator's ink dot (circle) changes completely and uniformly from red to yellow color. Compare the change to the Interpretation Guide printed on the chemical indicator and in this document.

5. If an unsatisfactory result is evidenced and the sterilization cycle fails, all materials contained in the load must be reviewed.

6. The chemical indicator is a legal document. Therefore, its result must be recorded in an appropriate book/control and be attached to the patient's medical record, according to the routine established in the institution and health legislation in force.

**CAUTION:** Do not use the **MAXITEST CHEMICAL INDICATOR - Type 4 (VH2O2)** to monitor sterilization cycles by saturated steam, ethylene oxide, dry heat, chemical vapors, radiation, or other sterilization processes that do not use plasma or hydrogen peroxide vaporization (VH2O2) as a sterilizing agent.

### VALIDITY/ANVISA REGISTRATION

The **MAXITEST CHEMICAL INDICATOR - Type 4 (VH2O2)** is valid for 02 years from the date of manufacture. Lot information, expiration date and other identification data are printed on the packaging and on each unit of the product.

Product exempt from registration with ANVISA, according to RDC 185/2001.

### STORAGE

Store away from light at temperatures between 15°C and 30°C and relative humidity between 30% and 60%. Do not freeze. Do not store in contact with floors and walls or near heat sources, sterilizing agents and other chemical products.

### BIO PROGRAM

Green consciousness is the logic that drives MAXIMED's business model. For this reason, the company is proud to have its people, facilities, processes, products and services incorporated into a sustainable system.

### GREEN PARTNER PROGRAM

To ensure that the components used in its products are harmless to the environment, MAXIMED incorporates the "Green Partner" into its processes, a program that approves the company's suppliers and partners according to appropriate environmental parameters and encourages the conduct of business that generates value over time.

### LEAD FREE

MAXIMED uses lead-free technologies and other heavy materials that dissolve when in contact with acidic water, contaminating water, air and soil, seriously harming human health.

### FREE FROM VOLATILE ORGANIC COMPOUNDS

MAXIMED uses technologies free from Volatile Organic Compounds that, when exposed to closed places, can worsen human health.

These are some of the actions in which MAXIMED reinforces its commitment to continuous improvement and the preservation of the environment. Copyright 2019. Clicktest, Maxitest and Maxipack are registered trademarks of Maximed Latin America. Reproduction in whole or in part without written permission is prohibited.