

MAXIPACK 3H

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PCD - PROCESS CHALLENGE DEVICE FLUORIMETRIC BIOLOGICAL INDICATOR 3H TYPE 5 CHEMICAL INTEGRATOR

PRODUCT DESCRIPTION

The **MAXXPACK CHALLENGE TEST 3H** was developed as an effective device to evaluate the performance of steam sterilization cycles and assist those responsible for MSC (Materials and Sterilization Center) in quality management in the reprocessing of medical materials and devices. System composed of a barrier that offers resistance to the action of the sterilizing agent associated with 01 **MAXXPACK Chemical Integrator - Type 5** that evaluates all the critical parameters of the cycle (steam quality, temperature and time), and 02 fast reading **CLICKTEST Fluorimetric Biological Indicators IBVF3H** in 03 hours, which makes it possible to prove microbial lethality.

Classified as a challenge test because it represents the most difficult load configuration to be sterilized. Indicated to monitor sterilization cycles from 121°C to 134°C **in loads that contain implantable health products** (ANVISA RDC 15:2012).

The **MAXXPACK CHALLENGE TEST 3H** is available in boxes of 10 units and complies with AAMI/ISO 11138-1, AAMI/ISO 11138-3, AAMI/ISO 11140-1, AAMI/ISO 17665-1 and AAMI ST 79 standards.

COMPOSITION

Package produced in impermeable paper with an external label containing a chemical process indicator (Type 1) in pink; an internal set of porous sheets that act as an important barrier system to the removal of air and the penetration of steam, complying with the standard established by the AAMI ST 79 regulation. It also has 01 **MAXXPACK Chemical Integrator - Type 5**, made of paper (medical grade) laminated on both sides, printed with steam thermochemical formula. Named as an integrator for carrying out a complete and integrated assessment of the critical parameters of the sterilization cycle (steam quality, time and temperature). For internal use and developed in accordance with the AAMI/ISO 11140-1 standard, it guarantees efficiency and safety in the application. Dimensions: 95mm (H) x 95mm (W). As a differential, the **MAXXPACK CHALLENGE TEST 3H** includes 02 **CLICKTEST IBVF3H STEAM Biological Indicators** made with: thermoplastic tube, stainless steel ACS, special filter, adhesive label with chemical process indicator (Type 1), minimum spore population of 1x10⁹ of Geobacillus stearothermophilus ATCC 7953 inoculated on a paper strip (spore-bearing), purple-colored culture medium contained in a wrapper and a thermoplastic lid with holes, added with a vapor permeable barrier. Complies with AAMI/ISO 11138 parts 1 and 3.

To endorse MAXXIMED's commitment to the conscientious use of resources and the preservation of the environment, **MAXXPACK CHALLENGE TEST 3H** is made entirely of recyclable materials and its compact size (115 mm X 95 mm) makes it easy to store. The unique, innovative design of the CLICKTEST Biological Indicator dispenses glass components and combined with an exclusive activation system, simplifies use and provides safe handling by preventing any risk of accidents.

INSTRUCTIONS FOR USE

1. Take a **MAXXPACK CHALLENGE TEST 3H** package and note the identification data: Auto-Reader Incubators number, cycle number and processing date.

2. Remove the **CLICKTEST IBVF3H** Biological Indicator (Positive Control) from the niche located at the bottom of the package and also identify it, noting on its label the Auto-Reader Incubators number, the cycle number and the processing date. **Reserve.**

3. Place the package and the materials to be sterilized in an area considered critical (of greatest challenge) of the Auto-Reader Incubators, usually located near the door and next to the drain. Make sure it is in the horizontal position and with the side that contains the Chemical Indicator (Type 1) in pink facing up.

4. Start the sterilization cycle corresponding to the load to be processed. After the end of the cycle, remove the **MAXXPACK CHALLENGE TEST 3H** package from the Auto-Reader Incubators, respecting the cooling time recommended by the manufacturer.

5. Confirm that the external chemical indicator (Type 1) contained in the package has changed from its pink (original) color to brown/black, confirming that it has been exposed to steam at 121°C or above.

Important: color change is not evidence that the sterilization cycle was successful. If there is no color change, it will be necessary to review the sterilization process.

6. Open the **MAXXPACK CHALLENGE TEST 3H** package through the perforation located on the back and remove the **MAXXPACK Chemical Integrator - Type 5** and the **CLICKTEST IBVF3H Fluorimetric Biological Indicator (TEST)** inside the package.

7. Verify that the type 5 Chemical Integrator has entirely changed from pink color to shades of brown to black. Compare the change to the Interpretation Guide printed on the strip and this Manual.

IMPORTANT: ACCEPT THE CHEMICAL INTEGRATOR WITH RESULTS EQUAL TO OR DARKER THAN THE REFERENCE (IN TONES FROM BROWN TO BLACK COLOR).

NOTE: Black color means the process was successful. However, the exposure time exceeded that required for efficient sterilization. These thermochemical ink shade variations occur due to the heating period of the Auto-Reader Incubators (heating ramp) until reaching the sterilization temperature. This heating time exposes the chemical integrators to sterilizing temperatures for a time much longer than the time calibrated in a resistometer (ISO 11140). Thus, variations in shades of brown and black are acceptable due to the non-uniformity of sterilization cycles and different load conditions.

8. If there is a failure in the process, the material must be processed again using a new **MAXXPACK CHALLENGE TEST 3H**.

9. Record the results according to current legislation and the institution's protocol.

10. Confirm that the external chemical indicator (Type 1), contained on the label of the **processed CLICKTEST Biological Indicator -IBVF3H** has changed from its pink (original) color to brown or black, confirming that there has been exposure to steam at 121°C or above.

Important: color change is not evidence that the sterilization cycle was successful. If there is no color change, it will be necessary to review the sterilization process.

11. Identify the **CLICKTEST Biological Indicator - IBVF3H STEAM** submitted to the sterilization process (processed), noting the Auto-Reader Incubators number, lot/load number and processing date.

12. Activate the 02 **CLICKTEST Biological Indicators - IBVF3H STEAM** (processed and positive control). To do this, hold the ampoule and press the lid down **until you hear a "click"**, breaking the retaining seal of the culture medium. Shake it **until all the liquid** (culture medium) **has completely descends** from the lid and comes into contact with the spore-bearing strip (see quick manual contained in this leaflet).

13. Ensure that the **fluorescence technology** incubator operates between 58°C and 62°C and that the programmed time is compatible with the biological indicators of 03 hours.

QUICK READ RESULT - FLUORIMETRIC METHODOLOGY

Incubate the 02 biological indicators (processed and positive control) for 03 hours. After the end of the incubation, observe:

a) if the incubator reader did **NOT** detect a fluorimetric reaction in the **processed indicator**, it characterizes a **negative result** and a successful sterilization cycle due to the inactivation of the spore population.

b) detecting fluorescence by the incubator means a **positive result** for the survival of microorganisms and a failed sterilization cycle.

14. Release the load/load for use only if the processed **Clicktest Biological Indicator - IBVF3H (STEAM)** presents a result "Negative - Absence of Microorganisms" and the control indicator results "Positive - Presence of microorganisms", confirming that the sterilization cycle was satisfactory and that the incubation was adequate.

Note: if the processed biological indicator shows a result "Positive - Presence of microorganisms" and, consequently, an irregularity in the sterilization cycle, the Hospital Infection Control Commission and Engineering/Maintenance must be notified to adopt the appropriate measures.

15. Detach the biological indicator labels (processed and positive control), attach them to the designated book/file and record the results obtained.

16. Dispose of the biological indicator tubes/ampoules following the sanitary regulations in your region. Positive indicators can be sterilized in cycles of 121°C for 15 minutes or 132°C for 10 minutes in gravitational Auto-Reader Incubators; or in cycles of 132°C for 04 minutes in vacuum pump-assisted Auto-Reader Incubators.

All the materials that make up the **MAXXPACK CHALLENGE TEST 3H**, including the biological indicators and the Type 5 chemical integrator are 100% recyclable, complying with the MAXXIMED institutional plan of ecological design.

NOTE: if necessary, it is possible to confirm the result obtained by rapid reading, complementing it with colorimetric technology (reading by changing the color of the culture medium) as follows:

RESULTS BY VISUAL READING

- **Do not perform steps 14, 15, and 16 above.**

- Check if the incubator is operating between 58°C and 62°C and incubate the 02 biological indicators for 24 hours to obtain the result by visual reading through colorimetric methodology (change of color of the culture medium). Since the presence of microorganisms can be detected before 24 hours, reading every 04 hours is recommended. Results:

a) maintenance of the purple color of the culture medium indicates the absence of viable microorganisms. Therefore, a successful sterilization process and **satisfactory result**.

b) change in the color of the culture medium, from purple to gray/greenish and later to yellow, indicates the presence/survival of microorganisms. Therefore, sterilization cycle failure and **unsatisfactory results**.

- Perform steps 14, 15 and 16 above

IMPORTANT: The 03 hours rapid fluorimetric reading and the optional 24-hour colorimetric reading are correlated with the 07 day incubation, according to the FDA's reduced incubation time protocol (RIT), and its results have a correlation of > or = to 97% of the cases.

CAUTION: Do not use the **MAXXPACK CHALLENGE TEST 3H** to monitor ethylene oxide sterilization cycles, dry heat, plasma/hydrogen peroxide vaporization, low temperature steam and formaldehyde, chemical vapors, radiation or other sterilization processes that do not use saturated steam as a sterilizing agent. Do not reuse biological indicators.

VALIDITY/ANVISA REGISTRATION

The **MAXXPACK CHALLENGE TEST 3H** package is valid for 02 years from the date of manufacture. Information about the batch, date of manufacture and expiration date are printed on the packaging of the **CLICKTEST IBVF3H Biological Indicators** and on the **MAXXPACK Chemical Integrator - Type 5**. Product exempt from ANVISA registration, according to RDC185/2001.

STORAGE

Store away from light at temperatures between 15°C and 30°C and relative humidity between 35% and 60%. Do not freeze. Do not store in contact with floors and walls or near heat sources, sterilizing agents, and other chemical products.